

Replacement Consumables for Millipore Lab Water Systems

Rephile offers a wide range of replacement cartridges and accessories for Millipore lab water systems.



with RFID tags

RephiDuo A Pack S2, replace Millipore PR0G0T0S2

Designed for water softening and to remove large particles, organic compounds and chlorine from tap water. The Cat. No. is RR800CP01.

Compatible with the following Millipore water systems:

- Elix Advantage
- Elix Reference
- Milli-Q Integral
- Elix Essential
- RiOs Essential
- AFS Essential



with RFID tags

RephiDuo A Pack, replace Millipore Q-Gard

Designed to remove ionic and organic contaminants down to trace levels from RO or distilled pure water.

Compatible with the following Millipore water systems:

- Milli-Q Advantage A10
- Milli-Q Reference

Repair services for Millipore PCB and power boards

Rephile provides repair services for failed control boards and power boards of various Millipore water systems and warrant full functionality of repaired boards for 6 months after service.

Millipore Cat. No.	Description	Millipore systems
ZF3000022	PCB ELIX/RiOs INTERFACE	RiOs (96 / Century) Elix (96 / Century)
ZF3000023	PCB MILLI-Q INTERFACE	Milli-Q (96 / Century)
ZF3000405 ZF3000487	MQ ADV_REF PWR BOARD	Milli-Q Advantage A10 Milli-Q Reference
ZF3000406	MQ ADV MAIN CONTROL BOARD	Milli-Q Integral (3/5/10/15) Milli-Q Direct (8/16) Milli-Q Advantage A10 Milli-Q Reference Elix Advantage (3/5/10/15) Elix Reference (3/5/10/15)
ZF3000433	MQ INT POWER SUPPLY BOARD	Milli-Q Integral (3/5/10/15) Milli-Q Direct (8/16) Elix Advantage (3/5/10/15) Elix Reference (3/5/10/15)

For any further information, please contact us by info@rephile.com.



A10 UV Lamp



Vent Filter



EDI



Final Filter

Cartridges

Millipore System	Millipore Cat. No.	RephiLe Part Number	Product Description
Milli-Q Academic Milli-Q Biocel Milli-Q Gradient Milli-Q Synthesis Milli-Q Element	QGARD00R1	RR400Q101	RephiDuo U Pack 1, replace QGARD00R1
	QGARD00D2	RR400Q201	RephiDuo U Pack 2, replace QGARD00D2
	QGARD00B1	RR400Q0B1	RephiDuo U Pack B1, replace QGARD00B1
	QTUM000IX	RRQTUM0IX	RephiSolo U Pack IX, replace QTUM000IX
	QTUM000EX	RRQTUM0EX	RephiSolo U Pack EX, replace QTUM000EX
	QTUMMPEIX	RRQTUMFIX	RephiSolo U Pack IX kit, replace QTUMMPEIX
	QTUMMPEEX	RRQTUMFEX	RephiSolo U Pack EX kit, replace QTUMMPEEX
Milli-Q Element	QGARD00B1	RR400Q0B1	RephiDuo U Pack B1, replace QGARD00B1
Milli-Q Element Q-POD Element	QTUM00ICP	RRQTUMICP	RephiSolo U Pack ICP, replace QTUM00ICP
Milli-Q Advantage A10 Milli-Q Reference	QGARDT1X1	RR800Q101	RephiDuo A Pack T1, replace QGARDT1X1
	QGARDT2X1	RR800Q201	RephiDuo A Pack T2, replace QGARDT2X1
Milli-Q Advantage A10 Milli-Q Integral Milli-Q Reference	QTUM0TIX1	M001IX	Replace QTUM0TIX1
	QTUM0TEX1	M01EX	Replace QTUM0TEX1
Milli-Q Direct	PR0G000T3	M001T3	Replace PR0G000T3
	QPAK00TIX	RR800QTIX	RephiDuo U Pack TIX, replace QPAK00TIX
	QPAK00TEX	RR800QTEX	RephiDuo U Pack TEX, replace QPAK00TEX
Milli-Q Plus Milli-Q PF Plus Milli-Q UF Plus Milli-Q UV Plus	CPMQ004R1	RR100Q101	RephiQuatro U Pack 1, replace CPMQ004R1
	CPMQK05R1	RR100Q1KT	RephiQuatro U Pack 1 kit, replace CPMQK05R1
	CPMQ004D2	RR100Q201	RephiQuatro U Pack 2, replace CPMQ004D2
	CPMQK05D2	RR100Q2KT	RephiQuatro U Pack 2 kit, replace CPMQK05D2
Elix Advantage Elix Reference Milli-Q Integral AFS Essential Elix Essential RiOs Essential	PR0G0T0S2	RR800CP01	RephiDuo A Pack S2, replace PR0G0T0S2
Elix 3/5/10 RiOs 5/8/16 Direct-Q 5	PR0G00002	RR400CP01	RephiDuo P Pack 2, replace PR0G00002
	PR0G000S2	RR400CPS2	RephiDuo P Pack S2, replace PR0G000S2
	PR0G00001	RR400CP00	RephiDuo P Pack 1, replace PR0G00001, 2/pk
Elix Clinical Elix Gulfstream Elix 20/35/70/100 RiOs 30/50/100/150/200	PR0GTL0S1	RR500CP01	RephiDuo L Pack TL1, replace PR0GTL0S1
	PR0GTLCS1	RR500CPC1	RephiDuo L Pack TL1, replace PR0GTL0S1
	PRPK0L0S1	RR500CPP1	RephiDuo L Pack L1, replace PRPK0L0S1
	QGARDTL01	RR500Q101	RephiDuo L Pack TL, replace QGARDTL01
Simplicity Simplicity 185 Synergy Synergy 185	SIMPAK0R1	RRSYNM0IX	RephiDuo S Pack 1, replace SIMPAK0R1
	SIMPAKKR1	RRSYNMFIX	RephiDuo S Pack 1 kit, replace SIMPAKKR1
	SIMPAK0D2	RRSYNM0DX	RephiDuo S Pack 2, replace SIMPAK0D2
	SIMPAKKD2	RRSYNMFDX	RephiDuo S Pack 2 kit, replace SIMPAKKD2
Simplicity (SMART) Simplicity UV	SIPK0SIA1	M001A1	Replace SIPK0SIA1
	SIPK0SIX2	M001X2	Replace SIPK0SIX2
Synergy (SMART) Synergy UV	SYPK0SIA1	M002A1	Replace SYPK0SIA1
	SYPK0SIX2	M002X2	Replace SYPK0SIX2
Direct-Q (SMART)	SPR00SIA1	RRSPR00A1	RephiDuo S Pack, replace SPR00SIA1
	SPR0LSIA1	RRSPR0LA1	RephiDuo S Pack, replace SPR0LSIA1
	SPR08SIA1	RRSPR08A1	RephiDuo S Pack, replace SPR08SIA1
RiOs (SMART)	SPR00S001	RRSPR0001	RephiSolo S Pack, replace SPR00S001
	SPR0LS001	RRSPR0L01	RephiSolo S Pack, replace SPR0LS001
	SPR08S001	RRSPR0801	RephiSolo S Pack, replace SPR08S001
RiOs DI (SMART)	SPR00SIA2	RRSPR00A2	RephiDuo S Pack, replace SPR00SIA2
AFS 10/10D/18/18D/30/30D AFS-25 EDI/AFS-50 EDI Milli-RO 3/5/8/12 Milli-RO 30/60/90 Milli-RX 20/45/75	CPR0P0402	RR100CP01	RephiQuatro P Pack, replace CPR0P0402
AFS 3/8/16	PR0G00003	RR400CP03	RephiSolo P Pack 3, replace PR0G00003
	CP1ALLRES	RR1ALLRES	RephiSolo Clinical Pack, replace CP1ALLRES
	CP2ALLRES	RR400Q401	RephiDuo Clinical Pack, replace CP2ALLRES
AFS 10/18/30 AFS 10D/18D/30D AFS-25 EDI/AFS-50 EDI	CP4ALLRES	RR100Q401	RephiQuatro Clinical Pack, replace CP4ALLRES
AFS Essential	CP2ALLREST	RR800Q501	RephiDuo Clinical Pack, replace CP2ALLREST
Milli-DI	CPDI000L1	RRDI000L1	RephiDuo DI Pack L, replace CPDI000L1
	PRPK00001	RRPRPK001	RephiDuo L Pack pretreatment, replace PRPK00001, 1/pk

Filters, UV Lamps and RO Membranes

Millipore System	Millipore Cat. No.	RephiLe Part Number	Product Description
Elix Milli-Q Integral	TANKMPK01	RATANKVN1	Tank Vent Filter with CO ₂ scavenger, replace TANKMPK01
RiOs Milli-Q Direct	TANKMPK02	RATANKVN2	Tank Vent Filter, replace TANKMPK02
S.D.S.	TANKVNT01	RATANKVT1	Vent Filter (S.D.S. Storage & Distribution System) , replace TANKVNT01
	TANKVNT02	RATANKVT2	Vent Filter with CO ₂ Trap (S.D.S. Storage System) , replace TANKVNT02
	ZFRE012F0	RATANKVF0	Vent filter 12" , replace ZFRE012F0
	ZFRE012FC	RATANKVFC	Vent filter 12" with CO ₂ trap, replace ZFRE012FC
Simplicity	SLFH02510	RJF324510	Vent Filter, replace SLFH02510, 10/pk
Ultrapure Water Systems POD	MPGP04001 MPGP02001 SIMFILTER	RAFFC0250	0.2 µm PES High Flus Capsule Filter
AFS Direct-Q 3/5 Elix (UV) / RiOs Elix Advantage / Elix Reference Elix Essential / RiOs Essential Milli-Q Direct Milli-Q Integral Milli-RO 3/5	CDRC00201	RAR000201	Small RO membrane, replace CDRC00201
	CDRC60201	RAR060201	RO Membrane for small RiOs& Elix, replace CDRC60201
	CDRC60202	RAR060202	RO Membranes for small RiOs& Elix, replace CDRC60202, 2/pk
	CDRC75202	RAR075202	RO membranes for Elix 15, replace CDRC75202, 2/pk
Milli-RO 6/8/10/12	CDRC012HG	RAR0012HG (Reuse existing housing)	RO membrane w/o housing, replace CDRC012HG
Elix 35/70/100 RiOs 50/100/150/200	CDRC050HG	RAR0050HG (Reuse existing housing)	RO membrane w/o housing, replace CDRC050HG
Milli-Q Gradient (Century) Milli-Q Synthesis (Century) Milli-Q Element Milli-Q Advantage A10 Milli-Q Direct Milli-Q Integral Milli-Q Reference	ZMQUVLP01	265021701	UV Lamp 185/254 nm, replace ZMQUVLP01
Milli-Q Gradient (96) Milli-Q Synthesis (96)	QUVLQSL01	RAUVLQSL1	Replacement UV Lamp, replace QUVLQSL01
Direct-Q UV Simplicity UV (SMART) Synergy UV (SMART)	SYN185UV1	265020601	UV Lamp 185 nm, replace SYN185UV1
Synergy UV (High Flow)	ZLUVWPM01	RAUV212B1	UV lamp 185/254 nm (>2012 model), replace ZLUVWPM01
Milli-Q (96/Century) A10 Milli-Q Advantage A10 Milli-Q Integral	ZFA10UVM1 ZFA10UV01	RAA10UVM1	UV Lamp for A10 TOC monitor
Elix UV Elix Advantage/Elix Refer- ence Elix Essential Milli-Q Direct Milli-Q Integral RiOs-DI 3 UV	ZLXUVLP01	265010601	UV Lamp 254 nm, replace ZLXUVLP01
Elix 20/35/70/100 RiOs 30/50/100/150/200 ASM (Cat. No.: TANKASMIN) AFS Essential	ZLXUVLPL1	265011001	UV Lamp 254 nm, replace ZLXUVLPL1
ASM (Cat. No.: TANKS50UV/ TANKS5LUV/ZFRES50UV/ ZFRES60UV)	ZFRES00UV	2650110C1	UV Lamp 254 nm, replace ZFRES00UV
Elix (UV) Elix Advantage Elix Reference Elix Essential Milli-Q Integral AFS 10E/15E	ZLX0EDI03	W3T101571	LPX03 CEDI Replacement Module, replace ZLX0EDI03
	ZLX0EDI05	W3T101572	LPX05 CEDI Replacement Module, replace ZLX0EDI05
	ZLX0EDI10	W3T101573	LPX10 CEDI Replacement Module, replace ZLX0EDI10
	ZLX0EDI15	W3T262701	LPX15 CEDI Replacement Module, replace ZLX0EDI15
30 PE reservoirs	FTPF05380	RAPP05380	30L Tank Level Sensor, replace FTFP05380
60 PE reservoirs	FTPF05381	RAPP05381	60L Tank Level Sensor, replace FTFP05381
100 PE reservoirs	FTPF06805	RAPP06805	100L Tank Level Sensor, replace FTFP06805
S.D.S. 200	FTPF07743	RAPP07743	200L Tank Level Sensor, replace FTFP07743
	ZF2000022	RAPP20022	200L Tank Level Sensor, replace ZF2000022

Pumps and Others

Millipore System	Millipore Cat. No.	Rephile Part Number	Product Description
Elix 3/5/10 RiOs 5/8/16	ZF3000000	RASP00000 (Reuse existing housing)	Small RiOs/Elix pump, 36VDC, replace ZF3000000
Milli-Q (96/Century)	ZF3000001	RASP00001 (Reuse existing housing)	Milli-Q 96/Century pump, 36VDC, replace ZF3000001
Elix Advantage Elix Reference Milli-Q Integral Milli-Q Direct	ZF3000431	RASP00431 (Reuse existing housing)	Integral Booster Pump, replace ZF3000431
Milli-Q Advantage Milli-Q Direct Milli-Q Integral Milli-Q Reference	ZF3000401	RASP00401 (Reuse existing housing)	Integral Recirc Pump, replace ZF3000401
Milli-Q (96/Century)	ZF3000004	RAUV00004	UV lamp ballast for MQ, replace ZF3000004
Elix UV	ZF3000131	RAUV00131	UV lamp ballast 6-10W, replace ZF3000131
Milli-Q Synergy	ZFMQ000PR	RAPR58561	Pressure Regulator
	ZLFI00001	RASDI0001	Rephile SDI Kit

Ultrapure Cartridges

- Filled with LeFil™ and OrganeFil™ media, which are composed of ion exchange resins and activated carbon
- Designed to remove specific types of water contaminants with large capacity

RephiSolo and RephiDuo U Pack



- Compatible with the following Millipore water systems:

- Milli-Q Academic
- Milli-Q Biocel
- Milli-Q Gradient
- Milli-Q Synthesis
- Milli-Q Element

RephiDuo S Pack



- Compatible with the following Millipore water systems:

- Synergy
- Simplicity

Pretreatment Filtration Cartridges

- Designed to remove most particles, organic compounds and chlorine from tap water to protect the RO membranes

RephiDuo P Pack

Contains activated carbon and depth filtration media in a polypropylene housing.



- Replace Millipore Progard
- Compatible with the following Millipore water systems:

- Elix 3/5/10
- RiOs 5/8/16
- Direct-Q 5

**Superior
Quality**

**Guaranteed
Performance**

RephiDuo L Pack

Suitable for large water systems and designed to ensure long service time.



- Replace Millipore Progard TL
- Compatible with the following Millipore water systems:

- Elix 20/35/70/100
- RiOs 30/50/100/150/200
- Elix Clinical
- Elix Gulfstream Clinical

**Highly
Cost Effective**

All rights reserved © Copyright 2023 RephiLe Bioscience, Ltd.

RephiLe Bioscience, Ltd.

Toll Free: +1-855-rephile (+1-855-737-4453)

Lit. No.: RFP0872210

Email: info@rephile.com